

HAAKE Viscometer 7 R+

Location of Machine: Composites Lab, RFM 1218

Location of SOP and Machine Operating & Safety Manual: Composites Lab website under resources; Composites Lab TRACS site; and Hardcopy near machine.

Emergency Contact:

- Call 911
- Call EHS & Risk Management at 512-245-3616
- Call Head Lab Technician, Dr. Ray Cook (office 512-245-2050)
- Call Dr. Jitendra S Tate (office 512-245-4872)

Before using this machine:

- You must have permission from Dr. Tate.
- You must have received formal training from technician or, trained research student (designated by Dr. Tate) related to machine safety and operation.
- You must read and understand **SOP and Machine Cleaning Manual.**
- You must use this machine under direct supervision of Dr. Tate or, Dr. Cook or, trained research student (designated by Dr. Tate).
- You must have signed “Lab Rules” document with Dr. Tate. This document must be signed every semester fall, spring, and summer (as applicable).
- If you do NOT follow above instructions you will be held responsible for your own safety and damages.

Safety Precautions:

Protective Equipment: Prior to performing this procedure, the following personal protective equipment must be obtained and ready for use: **Gloves, Safety Goggles, Face Mask, Lab Coat.**

Important Safeguards:

1. Prior to performing this procedure, the following safety equipment must be accessible and ready for use: (e.g. chemical fume hood, biological safety cabinet, laminar flow hood, chemical spill kits) **Fume hood**
2. All liquids should be drained to containers for chemical disposal and properly marked.
3. In the event that a hazardous material spill during this procedure, be prepared to clean with cleaner according to MSDS of materials used.



General Information

Viscosity is a measure of the resistance of a fluid which is being deformed by either shear stress or tensile stress. Viscosity and viscoelastic behavior of resins are very critical in composites manufacturing. Our lab is equipped with a classic rotational viscometer for the fast determination of viscosity as defined in ISO 2555 and more ASTM standards. For the determination of the rheological properties of a substance the measuring range can be changed by using a different speed or by a change of the spindle. The selection depends on the standard to be followed.

Specifications:

Operating Temperature Range: -20 to 100°C

Reproducibility: $\pm 1\%$




Hertz: 50-60Hz



Speed Range: 0.5-800rpm

Torque Range: 0.1-30mNm

Typical Testing Time: 1 min

Viscosity Range 1-1,000,000 mPas

Viscometer Operating procedure	
<ol style="list-style-type: none">1. Viscometer Assembly<ol style="list-style-type: none">a. Attach vertical pole to the base using the wrench provided in the caseb. Attach the viscometer to the vertical polec. Connect power cable to the viscometerd. <u>Optional (if needed):</u> Attach the thermometer to the back panele. <u>Optional (only for large beakers):</u> Attach the frame to the viscometer	
<ol style="list-style-type: none">2. Turn on the viscometer<ol style="list-style-type: none">a. Plug in the power cable to 110 V power outletb. Turn of switch on the back panel of the viscometer3. Calibration<ol style="list-style-type: none">a. After the viscometer is turned on the following message will appear: "Calibration is out of date, remove spindle and press enter"b. Make sure the spindle is not attachedc. Press enter	
<ol style="list-style-type: none">4. Starting up<ul style="list-style-type: none">• Selected required spindle type from the case (Usually: R2-largers)• Attach the spindle. <u>Note: Spindle has left-handed thread</u>• The machine is ready to use	
<ol style="list-style-type: none">5. Placing the beaker	

<ol style="list-style-type: none">a. Raise the viscometer to the highest level using a screw on the vertical poleb. Place beaker with liquid under the spindlec. Lower the viscometer until the spindle is submerged to <u>the spindle's mark</u>d. <u>Optional (if used)</u>: Submerge the thermometer to the liquid	
<p>6. Measuring</p> <ol style="list-style-type: none">a. Using “Enter”, “Up”, and “Down”: select spindle type and speedb. Press “Start”c. Use “Enter” to change the speedd. Use “Stop” to stop the process	
<p>7. <u>Notes</u></p> <ol style="list-style-type: none">a. Wait for 30-60 sec before readingb. Good results are in range 60%-80% of torquec. Viscometer with beep (and indicate “Error”) if under/over loadedd. Results depend on beaker and volume of liquide. Use the same beakers for comparison measurements	

8. Disassembly/Cleaning
 - a. Turn off the viscometer
 - b. Disassemble the viscometer
 - c. Clean the spindle and frame (if used) with cleaner according to MSDS of liquid used
 - d. Return the components to the case
 - e. Return the case to the drawer

