


Phosphorus 32



Characteristics

- Radioactive half-life: 14.3 days
- Decay mechanism: Beta emission
- Energy: $E_{\max} = 1709 \text{ keV}$
 $E_{\text{avg}} = 695 \text{ keV}$
- Contamination monitoring:
Thin window Geiger-Mueller detector
- Shielding: 1 cm lucite
- Dosimetry:
 - > micro-curie quantities – none,
 - > milli-curie quantities - Film badge, TLD ring, urinalysis

Decay Table

Days	0	1	2	3	4	5	6
7	1000	953	908	865	824	785	748
14	712	679	646	616	587	559	533
21	507	483	460	439	418	398	379
28	361	344	328	312	298	284	270
35	257	245	234	223	212	202	192
42	183	175	166	159	151	144	137
49	131	124	119	113	108	102	98
56	93	89	84	80	77	73	70
63	66	63	60	57	55	52	50

Dose Rate

- The dose rate on contact on the side of a 1 mCi delivery vial can be on the order of 1000 mrem/hr.
- If possible, avoid direct hand contact with vials and sources.
- When working with 100 μ Ci or more of P-32, work should be done behind a 1 cm lucite shield

Skin Dose

- One microcurie of P-32 in direct contact with 1 cm² of bare skin gives a dose rate to the skin of about 9 rem/hr.
- Always protect your skin and eyes when handling unsealed materials. Wear gloves, lab coats, safety glasses, and shoes

Detection

- A thin window G-M survey meter should always be available.
- A survey should be made immediately after use and any "hot spots" should be decontaminated.

Dosimetry

- Film badges and TLD rings should be worn for all P-32 work when handling 1 millicurie or more.

Waste

- Handle and store your radioactive waste carefully.
- The bottles for liquid waste should be placed in a secondary container (e.g. a bucket or tray) to contain spills or leaks.
- When more than a millicurie is involved, place 1 cm lucite in front of the container for shielding